

---

## Cortisol Mouse mAb

|               |  |
|---------------|--|
| Cat. No:      | V13401   |
| Name:         | Cortisol Mouse mAb   |
| Storage:      | Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. |
| Clone type:   | Monoclonal   |
| Clone No:     | 3A2  |
| Host/Source:  | Mouse  |
| Subtype:      | IgG2b  |
| Reactivity:   | Human  |
| Purity:       | >90% as determined by SDS-PAGE   |
| Format:       | Liquid   |
| Buffer:       | 0.01M PBS(pH7.4).  |
| Applications: | ELISA, ICA   |
| Antigen:      | Cortisol-KLH ( <a href="#">bs-25463P</a> )                                       |