

---

## Influenza A Nucleoprotein Rabbit mAb

Cat. No:	V22704
Name:	Influenza A Nucleoprotein Rabbit mAb
Storage:	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Monoclonal
Clone No:	2H7
Host/Source:	Rabbit
Subtype:	IgG
Reactivity:	(predicted:&nbsp;Influenza A virus)
Purity:	>90% as determined by SDS-PAGE
Format:	Liquid
Concentration:	1mg/ml
Buffer:	0.01M PBS(pH7.4)
Antigen:	Recombinant Influenza A Nucleoprotein
Pairing Information:	
V22703-V22704	