
Influenza A Nucleoprotein Mouse mAb

Cat. No:	V22708
Name:	Influenza A Nucleoprotein Mouse mAb
Storage:	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Monoclonal
Clone No:	8G2
Host/Source:	Mouse
Subtype:	IgG
Reactivity:	(predicted: Influenza A virus)
Purity:	>90% as determined by SDS-PAGE
Format:	Liquid
Buffer:	0.01M PBS(pH7.4)
Applications:	ELISA, ICA
Antigen:	Recombinant Influenza A Nucleoprotein
Pairing Information:	
V22708-V22712	