
Influenza A Nucleoprotein Mouse mAb

| | |
|---------------|--|
| Cat. No: | V22710 |
| Name: | Influenza A Nucleoprotein Mouse mAb |
| Storage: | Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. |
| Clone type: | Monoclonal |
| Clone No: | 4C2 |
| Host/Source: | Mouse |
| Subtype: | IgG |
| Reactivity: | (predicted: Influenza A virus) |
| Purity: | >90% as determined by SDS-PAGE |
| Format: | Liquid |
| Buffer: | 0.01M PBS(pH7.4) |
| Applications: | ELISA, ICA |
| Antigen: | Recombinant Influenza A Nucleoprotein |