
Influenza B Nucleoprotein Recombinant Rabbit mAb

Cat. No:	V22803
Name:	Influenza B Nucleoprotein Recombinant Rabbit mAb
Storage:	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Recombinant
Clone No:	2C7
Host/Source:	Rabbit
Subtype:	IgG
Reactivity:	Influenza B Nucleoprotein [B/Austria/1359417/2021 (Victoria)], Influenza B Nucleoprotein [B/Phuket/3073/2013 (Yamaggata)]
Purity:	>90% as determined by SDS-PAGE
Format:	Liquid
Buffer:	0.01M PBS(pH7.4)
Applications:	ELISA, ICA
Antigen:	Recombinant Influenza B Nucleoprotein: (bs-44026P)
Pairing Information:	
V22803-V22805	