

---

## Influenza B Nucleoprotein Mouse mAb

Cat. No:	V22806
Name:	Influenza B Nucleoprotein Mouse mAb
Clone type:	Monoclonal
Clone No:	10F4
Host/Source:	Mouse
Subtype:	IgG
Reactivity:	(predicted:&nbsp;Influenza B virus)
Purity:	>90% as determined by SDS-PAGE
Format:	Liquid
Buffer:	0.01M PBS(pH7.4)
Applications:	ELISA, ICA
Antigen:	Recombinant Influenza B Nucleoprotein
Pairing Information:	
V22807-V22806	