
Influenza B Nucleoprotein Mouse mAb

Cat. No:	V22807
Name:	Influenza B Nucleoprotein Mouse mAb
Storage:	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Monoclonal
Clone No:	8G9
Host/Source:	Mouse
Subtype:	IgG
Reactivity:	(predicted: Influenza B virus)
Purity:	>90% as determined by SDS-PAGE
Format:	Liquid
Buffer:	0.01M PBS(pH7.4)
Applications:	ELISA, ICA
Antigen:	Recombinant Influenza B Nucleoprotein
Pairing Information:	
V22807-V22810/V22806	