

www.bioss.com.cn sales@bioss.com.cn techsupport@bioss.com.cn 400-901-9800

## **AD7c-NTP Mouse mAb**

Cat. No: V3003

Name: AD7c-NTP Mouse mAb

Storage: Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Clone type: Monoclonal

Clone No: 16D12K Host/Source: Mouse Subtype: IgG2a

Reactivity: (predicted: Human)

Format: Liquid

Concentration: ≥1mg/ml

Buffer: 0.01M PBS (pH7.4) with 0.02% Proclin300.

Applications: ELISA

Antigen: KLH conjugated synthetic peptide derived from human Neuronal thread protein AD7c-

NTP (bs-41194P)