

---

## AD7c-NTP Mouse mAb

Cat. No:	V3003
Name:	AD7c-NTP Mouse mAb
Storage:	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Monoclonal
Clone No:	16D12K
Host/Source:	Mouse
Subtype:	IgG2a
Reactivity:	(predicted:&nbsp;Human)
Format:	Liquid
Concentration:	≥1mg/ml
Buffer:	0.01M PBS (pH7.4) with 0.02% Proclin300.
Applications:	ELISA
Antigen:	KLH conjugated synthetic peptide derived from human Neuronal thread protein AD7c-NTP ( <a href="#">bs-41187P</a> )