

## Recombinant human MMP2 Protein,C-His (HEK293)

Catalog Number: bs-47247P

Species: Human

AA Seq: 1-660/660

Tags: C-His

Activity: Not tested

Purity: > 85% as determined by SDS-PAGE

Purification: AC

Form: Liquid

Storage: 150 mM PBS

Stored at -70°C or -20°C for one year. Avoid repeated freeze/thaw cycles. Stored at 2 to 8 °C for 2 days after thawing.

### VALIDATION IMAGES



The purified product was run on 15% SDS-PAGE at reducing conditions. The gel was stained with Coomassie brilliant blue and destained. Note: The molecular weight of the protein as observed on SDS-PAGE under reducing conditions may differ from its predicted value, due to variations in electrophoretic mobilities.

### PRODUCT SPECIFIC PUBLICATIONS

**[IF=19]** Yuancai Ge. et al. Rapid and Scalable Preparation of Highly Uniform, Atomically Thin MSe<sub>2</sub> (M = Ti, Nb, Ta) Nanosheets as Ultra-Sensitive SERS Substrates for Lateral Flow Immunoassay. ADV FUNCT MATER. 2025 Jan;;2420786 LFIA ; . 10.1002/adfm.202420786