
Recombinant human SLC44A2 Protein,C-hFc & C-His (HEK293)

Catalog Number: bs-47262P

Species: Human

AA Seq: 61-230/706

Tags: C-hFc & C-His

Activity: Not tested

Purity: > 85% as determined by SDS-PAGE


Purification: AC

Form: Liquid

Storage: 150 mM PBS

Stored at -70°C or -20°C for one year. Avoid repeated freeze/thaw cycles. Stored at 2 to 8 °C for 2 days after thawing.

VALIDATION IMAGES



The purified product was run on 15% SDS-PAGE at reducing conditions. The gel was stained with Coomassie brilliant blue and destained. Note: The molecular weight of the protein as observed on SDS-PAGE under reducing conditions may differ from its predicted value, due to variations in electrophoretic mobilities.