

---

## Recombinant human WNT9A Protein, N-SUMO & N-His (HEK293)

Catalog Number: bs-47337P

Species: Human

AA Seq: 241-300/365

Tags: N-SUMO & N-His

Activity: Not tested

Purity: > 85% as determined by SDS-PAGE

Purification: AC


Form: Liquid

Storage: 150 mM PBS

Stored at -70°C or -20°C for one year. Avoid repeated freeze/thaw cycles. Stored at 2 to 8 °C for 2 days after thawing.

### VALIDATION IMAGES

---



The purified product was run on 15% SDS-PAGE at reducing conditions. The gel was stained with Coomassie brilliant blue and destained. Note: The molecular weight of the protein as observed on SDS-PAGE under reducing conditions may differ from its predicted value, due to variations in electrophoretic mobilities.