

Recombinant human ICAM1 Protein,C-His (HEK293)

Catalog Number: bs-47360P

Species: Human

AA Seq: 1-480/532

Tags: C-His

Activity: Not tested

Purity: > 85% as determined by SDS-PAGE


Purification: AC

Form: Liquid

Storage: 150 mM PBS

Stored at -70°C or -20°C for one year. Avoid repeated freeze/thaw cycles. Stored at 2 to 8 °C for 2 days after thawing.

VALIDATION IMAGES



Left: The purified protein was run on 15% SDS-PAGE at reducing conditions. The gel was then stained with Coomassie brilliant blue and destained. Right: Western Blotting using HRP conjugated 6xHis antibody. Note: The molecular weight of the protein as observed on SDS-PAGE under reducing conditions may differ from its predicted value, due to variations in electrophoretic mobilities.