

---

## **Recombinant human PLXNB1 Protein,C-hFc & C-His (HEK293)**

Catalog Number: bs-47383P

Species: Human

AA Seq: 691-840/2135

Tags: C-hFc & C-His

Activity: Not tested

Purity: > 85% as determined by SDS-PAGE

Purification: AC

Form: Liquid

Storage: 150 mM PBS

Stored at -70°C or -20°C for one year. Avoid repeated freeze/thaw cycles. Stored at 2 to 8 °C for 2 days after thawing.