

---

## **Recombinant human SCGB3A2 Protein,N-MBP & C-His (HEK293)**

Catalog Number: bs-47418P

Species: Human

AA Seq: 22-93/93

Tags: N-MBP & C-His

Activity: Not tested

Purity: > 85% as determined by SDS-PAGE

Purification: AC

Form: Liquid

Storage: 150 mM PBS

Stored at -70°C or -20°C for one year. Avoid repeated freeze/thaw cycles. Stored at 2 to 8 °C for 2 days after thawing.