
Recombinant human CCL26 Protein,N-SUMO & N-His (HEK293)

Catalog Number: bs-47453P

Species: Human

AA Seq: 31-90/94

Tags: N-SUMO & N-His

Activity: Not tested

Purity: > 85% as determined by SDS-PAGE


Purification: AC

Form: Liquid

Storage: 150 mM PBS

Stored at -70°C or -20°C for one year. Avoid repeated freeze/thaw cycles. Stored at 2 to 8 °C for 2 days after thawing.

VALIDATION IMAGES



Left: The purified protein was run on 15% SDS-PAGE at reducing conditions. The gel was then stained with Coomassie brilliant blue and destained.