

bs-4425R

[Primary Antibody]

PTD016 Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Rabbit	Isotype: IgG	Applications: WB (1:500-2000) ELISA (1:5000-10000)
Clonality: Polyclonal		
GeneID: 51136	SWISS: Q5M7Z0	Reactivity: (predicted: Human, Mouse, Rat, Horse)
Target: PTD016		
Immunogen: KLH conjugated synthetic peptide derived from human PTD016: 31-130/435.		
Purification: affinity purified by Protein A		Predicted MW.: 50 kDa
Concentration: 1mg/ml		Subcellular Location: Cytoplasm
Storage: Preservative: 0.02% Proclin300, Constituents: 1% BSA, 0.01M PBS, pH7.4. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.		
Background: RNFT1 (Ring Finger Protein, Transmembrane 1) is a Protein Coding gene. Diseases associated with RNFT1 include Sertoli Cell-Only Syndrome. An important paralog of this gene is RNFT2.		