

bs-10694R

[Primary Antibody]

Forchlorfenuron Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Rabbit

Isotype: IgG

Applications: ELISA (1:5000-10000)

Clonality: Polyclonal

Reactivity: Species independent

Target: Forchlorfenuron

Purification: affinity purified by Protein A

Concentration: 1mg/ml

Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.

Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Background: Forchlorfenuron is a plant growth regulator. It has been approved for use on kiwi fruit and grapes in the USA, and it has been associated with exploding watermelons in China.