

**bs-10694R**

**[ Primary Antibody ]**

## **Forchlorfenuron Rabbit pAb**



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

### **— DATASHEET —**

**Host:** Rabbit

**Isotype:** IgG

**Applications:** **ELISA** (1:5000-10000)

**Clonality:** Polyclonal

**Reactivity:** (predicted: Forchlorfenuron)

**Target:** Forchlorfenuron

**Purification:** affinity purified by Protein A

**Concentration:** 1mg/ml

**Storage:** 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.

Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

**Background:** Forchlorfenuron is a plant growth regulator. It has been approved for use on kiwi fruit and grapes in the USA, and it has been associated with exploding watermelons in China.