
Forchlorfenuron Rabbit pAb

Catalog Number: bs-10694R

Target Protein: Forchlorfenuron

Concentration: 1mg/ml

Form: Liquid

Host: Rabbit

Clonality: Polyclonal

Isotype: IgG

Applications: ELISA (1:5000-10000)

Reactivity: (predicted:Forchlorfenuron)

Purification: affinity purified by Protein A

Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.

Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Background: Forchlorfenuron is a plant growth regulator. It has been approved for use on kiwi fruit and grapes in the USA, and it has been associated with exploding watermelons in China.