

---

## Influenza B Nucleoprotein Mouse mAb

Cat. No:	V22814
Name:	Influenza B Nucleoprotein Mouse mAb
Storage:	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Monoclonal
Clone No:	8G2
Host/Source:	Mouse
Subtype:	IgG1
Reactivity:	Influenza B Nucleoprotein [B/Austria/1359417/2021 (Victoria)], Influenza B Nucleoprotein [B/Phuket/3073/2013 (Yamaggata)]
Purity:	>90% as determined by SDS-PAGE
Format:	Lyophilized or Liquid
Buffer:	0.01M PBS(pH7.4)
Applications:	ELISA, ICA
Antigen:	Recombinant Influenza B Nucleoprotein ( <a href="#">bs-44026P</a> )
Pairing Information:	
V22815-V22814	