

---

## Influenza B Nucleoprotein Mouse mAb

|                      |                                                                                                                           |
|----------------------|---------------------------------------------------------------------------------------------------------------------------|
| Cat. No:             | V22815                                                                                                                    |
| Name:                | Influenza B Nucleoprotein Mouse mAb                                                                                       |
| Storage:             | Shipped at 2-8°C.                                                                                                         |
| Clone type:          | Monoclonal                                                                                                                |
| Clone No:            | 5A1                                                                                                                       |
| Host/Source:         | Mouse                                                                                                                     |
| Subtype:             | IgM                                                                                                                       |
| Reactivity:          | Influenza B Nucleoprotein [B/Austria/1359417/2021 (Victoria)], Influenza B Nucleoprotein [B/Phuket/3073/2013 (Yamaggata)] |
| Purity:              | >90% as determined by SDS-PAGE                                                                                            |
| Format:              | Lyophilized or Liquid                                                                                                     |
| Buffer:              | 0.01M PBS(pH7.4)                                                                                                          |
| Applications:        | ELISA, ICA                                                                                                                |
| Antigen:             | Recombinant Influenza B Nucleoprotein ( <a href="#">bs-44026P</a> )                                                       |
| Pairing Information: | V22816&V22814-V22815; V22815-V22816                                                                                       |