
Influenza B Nucleoprotein Mouse mAb

Cat. No:	V22816
Name:	Influenza B Nucleoprotein Mouse mAb
Storage:	Shipped at 2-8°C.
Clone type:	Monoclonal
Clone No:	2F4
Host/Source:	Mouse
Subtype:	IgM
Reactivity:	Influenza B Nucleoprotein [B/Austria/1359417/2021 (Victoria)], Influenza B Nucleoprotein [B/Phuket/3073/2013 (Yamaggata)]
Purity:	>90% as determined by SDS-PAGE
Format:	Lyophilized or Liquid
Buffer:	0.01M PBS(pH7.4)
Antigen:	Recombinant Influenza B Nucleoprotein (bs-44026P)
Pairing Information:	V22815-V22814&V22816; V22816-V22814&V22815