
Influenza A Nucleoprotein Mouse mAb

Cat. No:	V22717
Name:	Influenza A Nucleoprotein Mouse mAb
Storage:	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Monoclonal
Clone No:	1B1
Host/Source:	Mouse
Subtype:	IgG2b
Reactivity:	Influenza A Nucleoprotein [A/Victoria/4897/2022 (H1N1)], Influenza A Nucleoprotein [A/Wisconsin/67/2022(H1N1)], Influenza A Nucleoprotein [A/Victoria/2570/2019(H1N1)], Influenza A Nucleoprotein [A/Wisconsin/588/2019(H1N1)], Influenza A Nucleoprotein [A/Darwin/9/2021(H3N2)], Influenza A Nucleoprotein [A/Darwin/6/2021(H3N2)], Influenza A Nucleoprotein [A/Massachusetts/18/2022(H3N2)], Influenza A Nucleoprotein [A/Thailand/1(KAN-1)/2004 (H5N1)], Influenza A Nucleoprotein [A/Hong_Kong/483/1997(H5N1)]
Purity:	>90% as determined by SDS-PAGE
Format:	Lyophilized or Liquid
Buffer:	0.01M PBS(pH7.4)
Applications:	ELISA, ICA
Antigen:	Recombinant Influenza A Nucleoprotein (bs-44023P)
Pairing Information:	V22719&V22713-V22717