
E2 Mouse mAb

| | |
|----------------------|--|
| Cat. No: | V13101 |
| Name: | E2 Mouse mAb |
| Storage: | Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. |
| Clone type: | Monoclonal |
| Clone No: | 2F3 |
| Host/Source: | Mouse |
| Subtype: | IgG1 |
| Reactivity: | Human |
| Purity: | >90% as determined by SDS-PAGE |
| Format: | Lyophilized or Liquid |
| Buffer: | 0.01M PBS(pH7.4) |
| Applications: | ELISA |
| Antigen: | E2-KLH (bs-4533PB) |
| Pairing Information: | |
| Capture- | bs-4533PB(ELISA) |
| Detection- | V13101(ELISA) |