
Testosterone Mouse mAb

Cat. No:	V13202
Name:	Testosterone Mouse mAb
Storage:	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Monoclonal
Clone No:	5D2
Host/Source:	Mouse
Subtype:	IgG1
Reactivity:	Human
Purity:	>90% as determined by SDS-PAGE
Format:	Lyophilized or Liquid
Buffer:	0.01M PBS(pH7.4)
Applications:	ELISA, ICA
Antigen:	Testosterone-KLH (bs-25802PB)
Pairing Information:	
Capture-bs-25802PB(ELISA)	
Detection-V13201& V13202 &V13203(ELISA);	
Capture-bs-25802PB(ICA)	
Detection-V13202 &V13203(ICA)	