
Progesterone Mouse mAb

Cat. No:	V13301
Name:	Progesterone Mouse mAb
Storage:	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Monoclonal
Clone No:	1C7
Host/Source:	Mouse
Subtype:	IgG1
Reactivity:	(predicted: Human)
Purity:	>90% as determined by SDS-PAGE
Format:	Lyophilized or Liquid
Buffer:	0.01M PBS(pH7.4)
Applications:	ELISA, ICA
Pairing Information:	
Capture-bs-25801PB(ELISA)	
Detection-V13301& V13302 (ELISA)	