

---

## Progesterone Mouse mAb

Cat. No:	V13302
Name:	Progesterone Mouse mAb
Storage:	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Monoclonal
Clone No:	2B8
Host/Source:	Mouse
Subtype:	IgG1
Reactivity:	Human
Purity:	>90% as determined by SDS-PAGE
Format:	Lyophilized or Liquid
Buffer:	0.01M PBS(pH7.4)
Applications:	ELISA, ICA
Antigen:	Progesterone-KLH ( <a href="#">bs-25801PB</a> )
Pairing Information:	
Capture	bs-25801PB(ELISA)
Detection	V13301& V13302 (ELISA);
Capture	bs-25801PB(ICA)
Detection	V13302 (ICA)