
PTPRB Rabbit pAb

Catalog Number: bs-4434R

Target Protein: PTPRB

Concentration: 1mg/ml

Form: Liquid

Host: Rabbit

Clonality: Polyclonal

Isotype: IgG

Applications: IHC-P (1:100-500), IHC-F (1:100-500), IF (1:100-500), ICC/IF (1:100-500), ELISA (1:5000-10000)

Reactivity: (predicted:Human, Mouse, Rat, Rabbit, Sheep, Cow, Dog)

Predicted MW: 222 kDa

Subcellular Cell membrane

Locations:

Entrez Gene: 5787

Source: KLH conjugated synthetic peptide derived from human PTPRB: 1551-1650/1997.

Purification: affinity purified by Protein A

Storage: Preservative: 0.02% Proclin300, Constituents: 1% BSA, 0.01M PBS, pH7.4.

Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.