
SCRN3 Rabbit pAb

Catalog Number: bs-4443R

Target Protein: SCRN3

Concentration: 1mg/ml

Form: Liquid

Host: Rabbit

Clonality: Polyclonal

Isotype: IgG

Applications: WB (1:500-2000), IHC-P (1:100-500), IHC-F (1:100-500), IF (1:100-500), ICC/IF (1:100-500),
ELISA (1:5000-10000)

Reactivity: (predicted:Human, Mouse, Rat, Pig, Sheep, Cow)

Predicted MW: 49 kDa

Entrez Gene: 79634

Swiss Prot: Q0VDG4

Purification: affinity purified by Protein A

Storage: Preservative: 0.02% Proclin300, Constituents: 1% BSA, 0.01M PBS, pH7.4.

Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Background: SCRN3 (Secernin 3) is a Protein Coding gene. Gene Ontology (GO) annotations related to this gene include dipeptidase activity. An important paralog of this gene is SCRN2.