

bs-4443R

[Primary Antibody]

SCRN3 Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

<p>Host: Rabbit</p> <p>Clonality: Polyclonal</p> <p>GeneID: 79634</p> <p>Target: SCRN3</p> <p>Purification: affinity purified by Protein A</p> <p>Concentration: 1mg/ml</p> <p>Storage: Preservative: 0.02% Proclin300, Constituents: 1% BSA, 0.01M PBS, pH7.4. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.</p> <p>Background: SCRN3 (Secernin 3) is a Protein Coding gene. Gene Ontology (GO) annotations related to this gene include dipeptidase activity. An important paralog of this gene is SCRN2.</p>	<p>Isotype: IgG</p> <p>SWISS: Q0VDG4</p>	<p>Applications: WB (1:500-2000) IHC-P (1:100-500) IHC-F (1:100-500) IF (1:100-500) ICC/IF (1:100-500) ELISA (1:5000-10000)</p> <p>Reactivity: (predicted: Human, Mouse, Rat, Pig, Sheep, Cow)</p> <p>Predicted MW.: 49 kDa</p> <p>Subcellular Location: Cytoplasm ,Nucleus</p>
---	--	---