

bs-19645R

[Primary Antibody]

OPA1L Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Rabbit	Isotype: IgG	Applications: WB (1:500-2000) ELISA (1:5000-10000)
Clonality: Polyclonal		
GeneID: 54838	SWISS: Q9NX94	Reactivity: (predicted: Human)
Target: OPA1L		
Immunogen: KLH conjugated synthetic peptide derived from human OPA1L: 201-300/342.		
Purification: affinity purified by Protein A		Predicted MW.: 38 kDa
Concentration: 1mg/ml		Subcellular Location: Cell membrane
Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.		

— SELECTED CITATIONS —

- **[IF=2.8]** Yunxiang Li. et al. Eif2s3y alleviated LPS-induced damage to mouse testis and maintained spermatogenesis by negatively regulating Adamts5. THERIOGENOLOGY. 2023 Aug;; IHC ;Mouse. 37586163