

---

## Rabbit Anti-Dog IgM, HRP conjugated

Catalog Number: bs-0373R-HRP

Target Protein: Rabbit Anti-Dog IgM

Concentration: 2.0 mg/ml

Form: Liquid

Host: Rabbit

Clonality: Polyclonal

Isotype: IgG

Applications: WB (1:1000-10000), IHC-P (1:100-500), IHC-F (1:100-1000), ELISA (1:1000-10000)

Reactivity: Dog

Purification: affinity purified by Protein A

Storage: 10 mM TBS (pH=7.4) with 1% BSA, 0.03% Proclin300 and 50% glycerol.

Store at -20°C for one year. Avoid repeated freeze/thaw cycles.