

**bs-0373R-HRP**

**[ Secondary Antibodies ]**



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

## **Rabbit Anti-Dog IgM, HRP conjugated**

### **— DATASHEET —**

<b>Host:</b> Rabbit	<b>Isotype:</b> IgG	<b>Applications:</b> <b>WB</b> (1:1000-10000)
<b>Clonality:</b> Polyclonal		<b>IHC-P</b> (1:100-500)
<b>Target:</b> Rabbit Anti-Dog IgM		<b>IHC-F</b> (1:100-1000)
<b>Purification:</b> affinity purified by Protein A		<b>ELISA</b> (1:1000-10000)
<b>Concentration:</b> 2.0 mg/ml		<b>Reactivity:</b> Dog
<b>Storage:</b> 10 mM TBS (pH=7.4) with 1% BSA, 0.03% Proclin300 and 50% glycerol. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.		