

**bs-0373R-Gold**

**[ Secondary Antibodies ]**



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

## **Rabbit Anti-Dog IgM, Gold conjugated**

### — DATASHEET —

**Host:** Rabbit

**Isotype:** IgG

**Clonality:** Polyclonal

**Target:** Rabbit Anti-Dog IgM

**Purification:** affinity purified by Protein A

**Concentration:** 0.4 mg/ml

**Storage:** 20mM TBS (pH8.0) with 1% BSA and 0.02% Proclin300.  
Store at 2-8°C for 3-6 months. Avoid repeated freeze/thaw cycles.

**Applications:** **ICA** (1:20-200)  
**IEM** (1:20-200)

**Reactivity:** Dog