
Rabbit Anti-Dog IgM, Gold conjugated

Catalog Number: bs-0373R-Gold

Target Protein: Rabbit Anti-Dog IgM

Concentration: 0.4 mg/ml

Form: Liquid

Host: Rabbit

Clonality: Polyclonal

Isotype: IgG

Reactivity: Dog

Purification: affinity purified by Protein A

Storage: 20mM TBS (pH8.0) with 1% BSA and 0.02% Proclin300.

Store at 2-8°C for 3-6 months. Avoid repeated freeze/thaw cycles.