
Rabbit Anti-Biotin, Cy5 conjugated

Catalog Number: bs-0311R-Cy5

Target Protein: Rabbit Anti-Biotin

Concentration: 2.0 mg/ml

Form: Liquid

Host: Rabbit

Clonality: Polyclonal

Isotype: IgG

Applications: IF (1:100-1000), Flow-Cyt (1:100-1000)

Excitation spectrum: 649nm

Emission spectrum: 670nm

Not yet tested in other applications.

Optimal working dilutions must be determined by the end user.

Purification: affinity purified by Protein A

Storage: 10 mM TBS (pH=7.4) with 1% BSA, 0.03% Proclin300 and 50% glycerol.

Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Background: Biotin is a water soluble vitamin, generally classified as a B complex vitamin, also called vitamin B4. After the initial discovery of biotin, nearly forty years of research were required to establish it as a vitamin. Biotin is required by all organisms but can only be synthesized by bacteria, yeasts, molds, algae, and some plant species. Biotin is required as prosthetic group of enzymes involved in incorporation of carbon dioxide into organic compounds. Biotin has a MW of 244 Da.