

**bs-3074R****[ Primary Antibody ]****phospho-MCSF Receptor (Tyr546) Rabbit pAb****BioSS**  
**ANTIBODIES**

www.bioss.com.cn


sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

**— DATASHEET —**

<b>Host:</b> Rabbit	<b>Isotype:</b> IgG	<b>Applications:</b> WB (1:500-2000)
<b>Clonality:</b> Polyclonal		<b>Reactivity:</b> Human (predicted: Mouse, Rat, Pig, Chicken, Dog, Horse)
<b>GeneID:</b> 1436	<b>SWISS:</b> P07333	<b>Predicted MW.:</b> 107 kDa
<b>Target:</b> MCSF Receptor (Tyr546)		<b>Subcellular Location:</b> Cell membrane
<b>Immunogen:</b> KLH conjugated Synthesised phosphopeptide derived from human CSF1R around the phosphorylation site of tyrosine 546: PK(p-Y)QV.		
<b>Purification:</b> affinity purified by Protein A		
<b>Concentration:</b> 1mg/ml		
<b>Storage:</b> 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.		
<b>Background:</b> This protein tyrosine kinase transmembrane receptor is the receptor for colony stimulating factor 1, a cytokine which controls the production, differentiation, and function of macrophages. This receptor mediates most if not all of the biological effects of this cytokine. Ligand binding activates the receptor kinase through a process of oligomerization and transphosphorylation. The encoded protein is a member of the CSF1/PDGF receptor family of tyrosine protein kinases and contains 5 immunoglobulin like C2 type domains. CD115 is expressed by cells of the monocytic lineage and by progenitor cells. Mutations in this gene have been associated with a predisposition to myeloid malignancy.		

**— VALIDATION IMAGES —**

Sample: Lane 1: Human U-87 MG cell lysates  
Lane 2: Human A431 cell lysates Primary: Anti-phospho-MCSF Receptor (Tyr546) (bs-3074R) at 1/300 dilution Secondary: IRDye800CW Goat Anti-Rabbit IgG at 1/20000 dilution Predicted band size: 107 kD Observed band size: 180 kD