bs-23294R

[Primary Antibody]

Bioss

Avi-Tag Rabbit pAb

www.bioss.com.cn sales@bioss.com.cn techsupport@bioss.com.cn 400-901-9800

- DATASHEET -

Host: Rabbit Isotype: IgG

Clonality: Polyclonal

Target: Avi-Tag

Purification: affinity purified by Protein A

Concentration: 1mg/ml

Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50%

Glycerol.

Shipped at 4°C. Store at -20°C for one year. Avoid repeated

freeze/thaw cycles.

Background: AviTag, a peptide allowing biotinylation by the enzyme BirA and so

the protein can be isolated by streptavidin (GLNDIFEAQKIEWHE)

Applications: WB (1:500-2000)

IHC-P (1:200-1000)
IHC-F (1:200-1000)
IF (1:200-1000)
Flow-Cyt (1ug/Test)
ICC/IF (1:100-500)
ELISA (1:5000-10000)

Reactivity: Avi Tag