


bsm-61414R**[Primary Antibody]****MIP3/CCL23 Recombinant Rabbit mAb****BioSS**
ANTIBODIES

www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —**Host:** Rabbit**Isotype:** IgG**Clonality:** Recombinant**GeneID:** 6368**SWISS:** P55773**Target:** MIP3/CCL23**Immunogen:** A synthesized peptide derived from human CCL23: 1-60.**Purification:** affinity purified by Protein A**Storage:** 10mM phosphate buffered saline(pH 7.4) with 150mM sodium chloride, 0.05% BSA, 0.02% Proclin300 and 50% glycerol. Store at 4°C for short term. Store at -20°C for long term. Avoid repeated freeze/thaw cycles.**Background:** Shows chemotactic activity for monocytes, resting T-lymphocytes, and neutrophils, but not for activated lymphocytes. Inhibits proliferation of myeloid progenitor cells in colony formation assays. This protein can bind heparin. Binds CCR1. CCL23(19-99), CCL23(22-99), CCL23(27-99), CCL23(30-99) are more potent chemoattractants than CCL23.**Applications:** **WB** (1:1000-2000)
ICC/IF (1:50-200)**Reactivity:** Human**Predicted
MW.:** 13**Subcellular
Location:** Secreted**— VALIDATION IMAGES —**

Western blot analysis of HepG2 cell lysate. Using MIP3/CCL23 (bsm-61414R) monoclonal antibody at 1:1000 dilution and 4°C overnight incubation. Followed by conjugated secondary antibody incubation at r.t. for 60 min.