

**bsm-61642R**

[ Primary Antibody ]



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

## phospho-Crkl (Tyr207) Recombinant Rabbit mAb

### — DATASHEET —

<b>Host:</b> Rabbit	<b>Isotype:</b> IgG	<b>Applications:</b> WB (1:1000-2000)
<b>Clonality:</b> Recombinant		<b>Reactivity:</b> Human
<b>GeneID:</b> 1399	<b>SWISS:</b> P46109	
<b>Target:</b> Crkl (Tyr207)		
<b>Immunogen:</b> A synthesized peptide derived from human CRKL around the phosphorylation site of Y207: AHA-pY-AQP.		<b>Predicted MW.:</b> 34
<b>Purification:</b> affinity purified by Protein A		<b>Subcellular Location:</b> Cytoplasm ,Nucleus
<b>Storage:</b> 10mM phosphate buffered saline(pH 7.4) with 150mM sodium chloride, 0.05% BSA, 0.02% Proclin300 and 50% glycerol. Store at 4°C for short term. Store at -20°C for long term. Avoid repeated freeze/thaw cycles.		
<b>Background:</b> May mediate the transduction of intracellular signals.		