

bsm-61783R**[Primary Antibody]****SMARCA5 Recombinant Rabbit mAb****BioSS**
ANTIBODIES

www.bioss.com.cn


sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Rabbit	Isotype: IgG	Applications: WB (1:1000-2000) ICC/IF (1:50-200)
Clonality: Recombinant		
GeneID: 8467	SWISS: O60264	Reactivity: Human, Mouse, Rat
Target: SMARCA5		
Immunogen: A synthesized peptide derived from human SNF2H: 1000-1052/1052.		
Purification: affinity purified by Protein A		Predicted MW.: 122
Storage: 10mM phosphate buffered saline(pH 7.4) with 150mM sodium chloride, 0.05% BSA, 0.02% Proclin300 and 50% glycerol. Store at 4°C for short term. Store at -20°C for long term. Avoid repeated freeze/thaw cycles.		Subcellular Location: Nucleus
Background: Helicase that possesses intrinsic ATP-dependent nucleosome-remodeling activity.		

— VALIDATION IMAGES —

Western blot analysis of MCF7 cell lysate. Using SMARCA5 (bsm-61783R) monoclonal antibody at 1:1000 dilution and 4°C overnight incubation. Followed by conjugated secondary antibody incubation at r.t. for 60 min.