
Rabbit Anti-Dog IgM, FITC conjugated

Catalog Number: bs-0373R-FITC

Target Protein: Rabbit Anti-Dog IgM

Concentration: 2.0 mg/ml

Form: Liquid

Host: Rabbit

Clonality: Polyclonal

Isotype: IgG

Applications: IF (1:100-1000), Flow-Cyt (1:100-1000)

Excitation spectrum: 495nm

Emission spectrum: 519nm

Not yet tested in other applications.

Optimal working dilutions must be determined by the end user.

Reactivity: Dog

Purification: affinity purified by Protein A

Storage: 10 mM TBS (pH=7.4) with 1% BSA, 0.03% Proclin300 and 50% glycerol.

Store at -20°C for one year. Avoid repeated freeze/thaw cycles.