

**bs-0373R-FITC**

**[ Secondary Antibodies ]**



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

## **Rabbit Anti-Dog IgM, FITC conjugated**

### **— DATASHEET —**

<b>Host:</b> Rabbit <b>Clonality:</b> Polyclonal <b>Target:</b> Rabbit Anti-Dog IgM <b>Purification:</b> affinity purified by Protein A <b>Concentration:</b> 2.0 mg/ml <b>Storage:</b> 10 mM TBS (pH=7.4) with 1% BSA, 0.03% Proclin300 and 50% glycerol. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.	<b>Isotype:</b> IgG	<b>Applications:</b> <b>IF</b> (1:100-1000) <b>Flow-Cyt</b> (1:100-1000) <b>ICC/IF</b> (1:100-1000) Excitation Spectrum: 495nm Emission spectrum: 519nm <b>Reactivity:</b> Dog
--	---------------------	---