

bs-0373R-Cy3

[Secondary Antibodies]



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

Rabbit Anti-Dog IgM, Cy3 conjugated

— DATASHEET —

Host: Rabbit	Isotype: IgG	Applications: IF (1:100-1000)
Clonality: Polyclonal		Flow-Cyt (1:100-1000)
Target: Rabbit Anti-Dog IgM		ICC/IF (1:100-1000)
Purification: affinity purified by Protein A		Excitation Spectrum: 552nm
Concentration: 2.0 mg/ml		Emission spectrum: 570nm
Storage: 10 mM TBS (pH=7.4) with 1% BSA, 0.03% Proclin300 and 50% glycerol. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.		Reactivity: Dog