


bsm-62392R**[Primary Antibody]****Fetuin A/AHSG Recombinant Rabbit mAb****Bioss**
ANTIBODIES

www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —**Host:** Rabbit**Isotype:** IgG**Clonality:** Recombinant**CloneNo.:** 14A13**GeneID:** 197**SWISS:** P02765**Target:** Fetuin A/AHSG**Immunogen:** A synthesized peptide derived from human Fetuin A: 11-57.**Purification:** affinity purified by Protein A**Storage:** 10mM phosphate buffered saline(pH 7.4) with 150mM sodium chloride, 0.05% BSA, 0.02% Proclin300 and 50% glycerol. Store at 4°C for short term. Store at -20°C for long term. Avoid repeated freeze/thaw cycles.**Background:** Promotes endocytosis, possesses opsonic properties and influences the mineral phase of bone. Shows affinity for calcium and barium ions.**Applications:** **WB** (1:500-2000)
IP (1:20-50)**Reactivity:** Human**Predicted
MW.:** 39**Subcellular
Location:** Secreted**— VALIDATION IMAGES —**

25 ug total protein per lane of various lysates
(see on figure) probed with Fetuin A/AHSG
monoclonal antibody, unconjugated
(bsm-62392R) at 1:1000 dilution and 4°C
overnight incubation. Followed by conjugated
secondary antibody incubation at r.t. for 60 min.