

bsm-62483R**[Primary Antibody]****Bioss**
ANTIBODIES

www.bioss.com.cn


sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

MRRF Recombinant Rabbit mAb**— DATASHEET —**

Host: Rabbit	Isotype: IgG	Applications: WB (1:1000-2000)
Clonality: Recombinant		IHC-P (1:100-200)
GeneID: 92399	SWISS: Q96E11	IHC-F (1:100-200)
Target: MRRF		IF (1:100-200)
Immunogen: A synthesized peptide derived from human MRRF: 235-262.		Reactivity: Human, Mouse
Purification: affinity purified by Protein A		
Storage: 10mM phosphate buffered saline(pH 7.4) with 150mM sodium chloride, 0.05% BSA, 0.02% Proclin300 and 50% glycerol. Store at 4°C for short term. Store at -20°C for long term. Avoid repeated freeze/thaw cycles.		Predicted MW.: 29
Background: Responsible for the disassembly of ribosomes from messenger RNA at the termination of mitochondrial protein biosynthesis. Acts in collaboration with GFM2. Promotes mitochondrial ribosome recycling by dissolution of intersubunit contacts.		Subcellular Location: Cytoplasm

— VALIDATION IMAGES —

Western blot analysis of HeLa cell lysate. Using MRRF (bsm-62483R) monoclonal antibody at 1:1000 dilution and 4°C overnight incubation. Followed by conjugated secondary antibody incubation at r.t. for 60 min.