

bsm-62494R

[Primary Antibody]

ACAT2 Recombinant Rabbit mAb

Bioss
ANTIBODIES

www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Rabbit

Isotype: IgG

Clonality: Recombinant

GeneID: 39

SWISS: Q9BWD1

Target: ACAT2

Immunogen: A synthesized peptide derived from human ACAT2: 200-397.

Purification: affinity purified by Protein A

Storage: 10mM phosphate buffered saline(pH 7.4) with 150mM sodium chloride, 0.05% BSA, 0.02% Proclin300 and 50% glycerol.
Store at 4°C for short term. Store at -20°C for long term. Avoid repeated freeze/thaw cycles.


Background: Involved in the biosynthetic pathway of cholesterol.

Applications: WB (1:1000-2000)

Reactivity: Human, Mouse, Rat

Predicted MW.: 41

— VALIDATION IMAGES —



Western blot analysis of K562 cell lysate. Using ACAT2 (bsm-62494R) monoclonal antibody at 1:1000 dilution and 4°C overnight incubation. Followed by conjugated secondary antibody incubation at r.t. for 60 min.