

bsm-62512R**[Primary Antibody]****BioSS**
ANTIBODIES

www.bioss.com.cn


sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

NCR2 Recombinant Rabbit mAb**— DATASHEET —**

Host: Rabbit	Isotype: IgG	Applications: WB (1:500-2000) ICC/IF (1:50-200)
Clonality: Recombinant		
GeneID: 9436	SWISS: O95944	Reactivity: Human
Target: NCR2		
Immunogen: A synthesized peptide derived from human NK p44: 25-60.		
Purification: affinity purified by Protein A		Predicted MW.: 31
Storage: 10mM phosphate buffered saline(pH 7.4) with 150mM sodium chloride, 0.05% BSA, 0.02% Proclin300 and 50% glycerol. Store at 4°C for short term. Store at -20°C for long term. Avoid repeated freeze/thaw cycles.		Subcellular Location: Cell membrane
Background: Cytotoxicity-activating receptor that may contribute to the increased efficiency of activated natural killer (NK) cells to mediate tumor cell lysis.		

— VALIDATION IMAGES —

Western blot analysis of K562 cell lysate. Using NCR2 (bsm-62512R) monoclonal antibody at 1:1000 dilution and 4°C overnight incubation. Followed by conjugated secondary antibody incubation at r.t. for 60 min.