

bsm-62590R**[Primary Antibody]****ENDOGL1 Recombinant Rabbit mAb****BioSS**
ANTIBODIES

www.bioss.com.cn


sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Rabbit	Isotype: IgG	Applications: WB (1:1000-2000) IP (1:20-50)
Clonality: Recombinant		
GeneID: 9941	SWISS: Q9Y2C4	Reactivity: Human
Target: ENDOGL1		
Immunogen: A synthesized peptide derived from human EXOG: 350-368.		
Purification: affinity purified by Protein A		
Storage: 10mM phosphate buffered saline(pH 7.4) with 150mM sodium chloride, 0.05% BSA, 0.02% Proclin300 and 50% glycerol. Store at 4°C for short term. Store at -20°C for long term. Avoid repeated freeze/thaw cycles.		
Background: Endo/exonuclease with nicking activity towards supercoiled DNA, a preference for single-stranded DNA and 5'-3' exonuclease activity.		
		Predicted MW.: 41
		Subcellular Location: Cell membrane

— VALIDATION IMAGES —

Western blot analysis of HepG2 cell lysate. Using ENDOGL1 (bsm-62590R) monoclonal antibody at 1:1000 dilution and 4°C overnight incubation. Followed by conjugated secondary antibody incubation at r.t. for 60 min.